

## Thermophilic Reverse Transcriptase

### Product Name and Catalog number

Thermophilic Reverse Transcriptase:  
 Taq-Standard-SYBR Green, Cat. # W140-TS  
 Taq-Fast-SYBR Green, Cat. # W140-TF  
 Taq-Probe-SYBR Green, Cat. # W140-TP  
 Size: 200U for 200 reactions (No ROX)

### Intended Use

This reverse transcriptase or RNA-dependent DNA polymerase is used for real time-fluorescence RT-PCR with either TaqMan probe or SYBR Green dye. Research use only.

**Please Choose a DNA polymerase and associated 5x RT-PCR Buffer before you start:**

Depending on used DNA polymerases, the kit contains different buffers.

**Table 1. Different DNA polymerases with related buffers**

DNA polymerase	Fluorescence available	5x Buffer
<i>Taq</i> polymerase <sup>a</sup>	SYBR Green	<b>A</b>
<i>Taq-Fast</i> polymerase <sup>b</sup>	SYBR Green	<b>B</b>
<i>Taq-Probe</i> <sup>c</sup> polymerase	TaqMan Probe	<b>C</b>

### Footnotes of Table 1

<sup>a</sup> *Taq* DNA polymerase is not included in the kit.

<sup>b</sup> *Taq-Fast* polymerase, Cat.# W148, can extend more than 300 bases with short PCR cycling program.

<sup>c</sup> *Taq-Probe* polymerase, Cat.# W145, a specially engineered enzyme for S-shaped curve, is not included.

Part I is for *Taq* DNA polymerase with Taq-Standard SYBR Green Dye, W140-TS, page 2-3.

Part II is for *Taq-Fast* polymerase with SYBR Green Dye, W140-TF, page 4-5.

Part III is for *Taq-Probe* polymerase with TaqMan probe, W140-TP, page 6-7.

## Part I: Thermophilic Reverse Transcriptase

### Thermophilic Reverse Transcriptase for SYBR Green Dye together with Taq DNA polymerase:

Cat.# **W140-TS**, 200U for 200 reactions x 20µl Taq-Standard-SYBR Green (No ROX)

#### Intended Use

This reverse transcriptase or RNA-dependent DNA polymerase is used for real-time-fluorescent RT-PCR with SYBR Green dye.

#### Features

- The reverse transcriptase is a thermophilic type A polymerase (US patent pending) with optimal temperatures at 60-62°C (Tables 1-2).
- The RTase is easily heat-inactivated at ≥90°C for 1min.
- The enzyme efficiently synthesizes a complementary DNA strand on RNA template from a gene-specific primer, ≤1 unit per 20µL of singleplex reaction.
- The enzyme reversely-transcribes single digit copies of target RNA molecules consistently.
- The concentrations of the primers are variable depending on assay designs and thermocycling protocols (Table 4).
- The preferable PCR product size is ≤150bp.
- It is particularly suitable for one-step real-time qualitative and quantitative RT-PCR.

**Table 1. Enzyme properties**

Reverse transcriptase activity	Yes
5'-3' DNA polymerase activity	Yes
5'-3' exonuclease activity	Yes
3'-5' exonuclease activity	No
Incorporation of modified nucleotides	Yes, e.g., dUTP, fluorescence dye-labeled dNTPs

**Table 2. Optimal parameters**

Optimal RTase amount	0.5-1U/20µl reaction
Optimal temperature	60-62°C
Optimal incubation time	5-10min
Heat-Inactivated	≥90°C for 2min
Mg <sub>2</sub> concentration	1.5mM

dNTP concentration	Each 200uM
Primer concentration	Each ≥0.15µM, depending on primer design and thermocycling
Product size	Preferred ≤150bp

#### Unit Definition

One unit of the enzyme catalyzes incorporation of 1nmol of deoxyribonucleotides into polynucleotide in 10min.

#### Production Source

*E. coli* strain

**Table 3. Kit Contents**

Component	Amount
RTase, 10U/µl	200U/500U/ 2500U
5x Buffer A	1mL/2x1mL/10mL
User manual	1

#### Ship and Store

The kit can be transported below 4°C for up to 3 days. It can be stored at -20°C for 24 months.

#### Setup Reaction and Thermocycling

1. Thaw 5x Buffer and other reaction components at room temperature, mix each component, centrifuge, and then place on ice.
2. Determine the total volume for the number of reactions, add 5-10% extra volume, and prepare assay mix of all components except RNA template. Mix the assay mix, centrifuge, and then place on ice.
3. Aliquot the assay mix into PCR tubes or plates.
4. Add RNA template to PCR tubes or plate.
5. Seal tubes with flat, optically transparent caps or seal plates with optically transparent film.
6. Mix and then briefly centrifuge the tubes or plate.
7. Program PCR instrument with indicated thermo-cycling protocol.
8. Load PCR tubes or plates and start to run.
9. Perform data analysis according to the PCR instrument instructions.

**Table 4. Set up a 20µl of RT-PCR reaction.**

Component	Amount or final concentration
5x Buffer A	4ul
dNTPs	Each 200uM
Target specific primers <sup>a</sup>	Each 0.15-0.2uM

RTase, 10U/μl	0.5-1U
<b>Taq polymerase<sup>c</sup></b>	1U
RNA template <sup>d</sup>	As low as single digit copies of target RNA
Nuclease-free H <sub>2</sub> O	To a total volume of 20μl

**Footnotes of Table 4**

<sup>a</sup> The primer T<sub>m</sub> should be designed ≥60°C, preferably between 62°C to 65°C, using primer3 software for high efficiency and specificity.

<sup>d</sup> RNA templates should be extracted by a qualified silica-based kit and eluted with low EDTA TE buffer (10mM Tris-HCl, 0.1mM EDTA, pH 8.0-8.3).

**Table 5. Compatible instruments**

RT-PCR Instrument	ROX required by instrument	Passive dye setup
Bio-Rad® iQ™5, CFX96, CFX384, Opticon Roche Lightcycler® Qiagen Rotor-Gene™ Eppendorf Mastercycler® Cepheid® SmartCycler®	Not recommended	Not necessary
Applied Biosystems® 7500, 7500 Fast, QuantStudio™, ViiA7™, Agilent Mx™	Low ROX (50nM final concentration)	Turn on ROX passive reference dye button
Applied Biosystems® 5700, 7000, 7300, 7700, 7900, 7900HT, 7900HT Fast, StepOne™, StepOnePlus™	High ROX (500nM final concentration)	Turn on ROX passive reference dye button

**Table 6. Standard thermocycling protocol**

Stage	Temperature	Period	Number of cycles
I	60°C	10min	1
II	95°C	2min	1
III	95°C	10sec	35-40
	60°C, signal acquisition	60sec	

**Table 7. Three-Step Thermocycling Protocol**

Stage	Temperature	Period	Number of cycles
I	60°C	10min	1
II	95°C	2min	1
III	95°C	10sec	35-40
	60°C	30sec	
	68-72°C, signal acquisition	30sec	

**Footnotes of Tables 6 and 7**

The three-step thermocycling protocol in Table 7 increases overall DNA polymerase activity by 50%, a more effective protocol than Table 6.

The primer concentration used in Tables 6 and 7 is typically 0.15-0.2uM.

## Part II: Thermophilic Reverse Transcriptase

### Product Name and Catalog Number

Thermophilic Reverse Transcriptase for **SYBR Green Dye together with Taq-Fast polymerase**

Cat.#: **W140-TF**, 200U for 200 reactions x 20µL (No ROX)

### Intended Use

This reverse transcriptase or RNA-dependent DNA polymerase is used for real time-fluorescence RT-PCR with SYBR Green dye.

### Features

- The reverse transcriptase is a thermophilic type A polymerase (US patent pending) with optimal temperatures at 60-62°C (Tables 1-2).
- The RTase is easily heat-inactivated at ≥90°C for 1min.
- The enzyme efficiently synthesizes a complementary DNA strand on RNA template from a gene-specific primer, ≤1 unit per 20µL of singleplex reaction.
- The enzyme reversely-transcribes single digit copies of target RNA molecules consistently.
- The concentrations of the primers are variable depending on assay designs and thermocycling protocols (Table 4).
- The preferred PCR product size is ≤150bp.
- It is particularly suitable for one-step real-time qualitative and quantitative RT-PCR.

**Table 1. Enzyme properties**

Reverse transcriptase activity	Yes
5'-3' DNA polymerase activity	Yes
5'-3' exonuclease activity	Yes
3'-5' exonuclease activity	No
Incorporation of modified nucleotides	Yes, e.g., dUTP, fluorescence dye-labeled dNTPs

**Table 2. Optimal parameters**

Optimal RTase amount	≤1U/20µl reaction
Optimal temperature	60-62°C
Optimal incubation time	5-10min

Heat-Inactivated	≥90°C for 1min
Mg <sub>2</sub> concentration	1.5mM
dNTP concentration	Each 200uM
Primer concentration	Each ≥0.15µM, depending on primer design and thermocycling
Product size	Preferred ≤150bp

### Unit Definition

One unit of the enzyme catalyzes incorporation of 1nmol of deoxyribonucleotides into polynucleotide in 10min.

### Production Source

*E. coli* strain

**Table 3. Kit Contents**

Component	Amount
RTase, 10U/µl	200U/500U, 2500U
5x Buffer <b>B</b> (SYBR Green)	1mL/2x1mL/10mL
User manual	1

**Notes: Taq-Fast polymerase, Cat.# W148**, which can extend more than 300 bases with short cycling program, and it is not included in the kit.

### Transportation and Storage

The kit can be transported below 4°C for up to 3 days. It can be stored at -20°C for 24 months.

### Setup Reaction and Thermocycling

1. Thaw 5x Buffer and other reaction components at room temperature, mix each component, centrifuge, and then place on ice.
2. Determine the total volume for the number of reactions, add 5-10% extra volume, and prepare assay mix of all components except RNA template. Mix the assay mix, centrifuge and then place on ice.
3. Aliquot the assay mix into PCR tubes or plates.
4. Add RNA template to PCR tubes or plate.
5. Seal tubes with flat, optically transparent caps or seal plates with optically transparent film.
6. Mix and then briefly centrifuge the tubes or plate.
7. Program PCR instrument with indicated thermo-cycling protocol.
8. Load PCR tubes or plates and start to run.
9. Perform data analysis according to the PCR instrument instructions.

**Table 4. Set up a 20µl of RT-PCR reaction.**

Component	Amount or final concentration
5x Buffer <b>B</b>	4ul
dNTPs	Each 200uM
Target specific primers <sup>a</sup>	Each ≥0.15uM
RTase, 10U/µl	≤1U
<b>Taq-Fast polymerase<sup>b</sup></b>	1U
RNA template <sup>c</sup>	As low as single digit copies of target RNA
Nuclease-free H <sub>2</sub> O	To a total volume of 20µl

**Footnotes of Table 4**

<sup>a</sup> The primer T<sub>m</sub> should be designed ≥60°C, preferably between 62°C to 65°C, using primer3 software for high efficiency and specificity.

<sup>b</sup> Taq-Fast polymerase, which can extend more than 300 bases with a short cycling program, is not included in the kit.

<sup>c</sup> RNA templates should be extracted by a qualified silica-based kit and eluted with low EDTA TE buffer (10mM Tris-HCl, 0.1mM EDTA, pH 8.0-8.3).

**Table 5. Compatible instruments**

RT-PCR Instrument	ROX required by instrument	Passive dye setup
Bio-Rad® iQ™ 5, CFX96, CFX384, Opticon Roche Lightcycler® Qiagen Rotor-Gene™ Eppendorf Mastercycler® Cepheid® SmartCycler®	Not recommended	Not necessary
Applied Biosystems® 7500, 7500 Fast, QuantStudio™, ViiA7™, Agilent Mx™	Low ROX (50nM final concentration)	Turn on ROX passive reference dye button
Applied Biosystems® 5700, 7000, 7300, 7700, 7900, 7900HT, 7900HT Fast, StepOne™, StepOnePlus™	High ROX (500nM final concentration)	Turn on ROX passive reference dye button

**Table 6. Standard thermocycling protocol**

Stage	Temperature	Period	Number of cycles
I	60°C	10min	1
II	95°C	2min	1
III	95°C	10sec	35-40
	60°C, signal acquisition	60sec	
IV	60°C to 95°C	Various	1

**Table 7. Three-Step Thermocycling Protocol**

Stage	Temperature	Period	Number of cycles
I	60°C	10min	1
II	95°C	2min	1
III	95°C	10sec	35-40
	60°C	30sec	
	68-72°C, signal acquisition	30sec	
IV	68-72°C to 95°C	Various	1

**Footnotes of Tables 6 and 7**

The three-step thermocycling protocol in Table 7 increases overall DNA polymerase activity by 50%, a more effective protocol than Table 6.

The primer concentration used in Tables 6 and 7 is typically 0.15-0.2uM.

**Table 8. Fast thermocycling protocol**

Stage	Temperature	Period	Number of cycles
I	60°C	5-10min	1
II	95°C	1min	1
III	95°C	5sec	35-40
	60°C, signal acquisition	30sec	
IV	60°C to 95°C	Various	1

**Footnotes of Table 8**

The product size for the fast two-step thermocycling protocol is preferred to be less than 90bp.

The primer concentration used is typically between 0.4uM and 0.9uM.

## Thermophilic Reverse Transcriptase

### Product Name and Catalog number

Thermophilic Reverse Transcriptase for **TaqMan probe with Taq-Probe polymerase**

Cat. #: **W140-TP**, 200U for 200 reactions x 20 $\mu$ L (No ROX)

### Intended Use

This reverse transcriptase or RNA-dependent DNA polymerase is used for real-time-fluorescent RT-PCR with TaqMan probe.

### Features

- The reverse transcriptase is a thermophilic type A polymerase (US patent pending) with optimal temperatures at 60-62°C (Tables 1-2).
- The RTase is easily heat-inactivated at  $\geq 90^\circ\text{C}$  for 1min.
- The enzyme efficiently synthesizes a complementary DNA strand on RNA template from a gene-specific primer,  $\leq 1$  unit per 20 $\mu$ L of singleplex reaction.
- The enzyme reversely-transcribes single digit copies of target RNA molecules consistently.
- The concentrations of the primers are variable depending on assay designs and thermocycling protocols (Table 4).
- The preferred PCR product size is  $\leq 150\text{bp}$ .
- It is particularly suitable for one-step real-time qualitative and quantitative RT-PCR.

**Table 1. Enzyme properties**

Reverse transcriptase activity	Yes
5'-3' DNA polymerase activity	Yes
5'-3' exonuclease activity	Yes
3'-5' exonuclease activity	No
Incorporation of modified nucleotides	Yes, e.g., dUTP, fluorescence dye-labeled dNTPs

**Table 2. Optimal parameters**

Optimal RTase amount	Singleplex: 1 pair of primers: 0.5-1U/20 $\mu$ L reaction
	Multiplex: 2-4 pairs of primers: 1-2U/20 $\mu$ L reaction
Optimal temperature	60-62°C

Optimal incubation time	5-10min
Heat-Inactivated	$\geq 90^\circ\text{C}$ for 2min
Mg <sub>2</sub> concentration	1.5-3mM
dNTP concentration	Each 200uM
Primer concentration	Each $\geq 0.15\mu\text{M}$ , depending on primer design and thermocycling
Probe concentration	Each 0.15-0.25uM
Product size	Preferred $\leq 150\text{bp}$

### Unit Definition

One unit of the enzyme catalyzes incorporation of 1nmol of deoxyribonucleotides into polynucleotide in 10min.

### Production Source

*E. coli* strain

**Table 3. Kit Contents**

Component	Amount
RTase, 10U/ $\mu$ L	200U/500U/ 2500U
5x Buffer C	1mL/2x1mL/10mL
User manual	1

**Notes:** Taq-Probe polymerase, Cat. # W145, a specially engineered enzyme for S-shaped curve, is not included.

### Transportation and Storage

The kit can be transported below  $4^\circ\text{C}$  for up to 3 days. It can be stored at  $-20^\circ\text{C}$  for 24 months.

### Setup Reaction and Thermocycling

- Thaw 5x Buffer and other reaction components at room temperature, mix each component, centrifuge and then place on ice.
- Determine the total volume for the number of reactions, add 5-10% extra volume, and prepare assay mix of all components except RNA template. Mix the assay mix, centrifuge and then place on ice.
- Aliquot the assay mix into PCR tubes or plate.
- Add RNA template to PCR tubes or plate.
- Seal tubes with flat, optically transparent caps or seal plates with optically transparent film.
- Mix and then briefly centrifuge the tubes or plate.
- Program PCR instrument with indicated thermo-cycling protocol.
- Load PCR tubes or plates and start to run.
- Perform data analysis according to the PCR instrument instructions.



**Table 4. Set up a 20µl of RT-PCR reaction**

Component	Amount or final concentration
5x Buffer C	4ul
dNTPs	Each 200uM
Target specific primers <sup>a</sup>	Each 0.15-0.2uM
Probes <sup>b</sup>	Each 0.15-0.25uM
RTase, 10U/µl	Singleplex: 0.5-1U
	Multiplex: 1-2U
<b>Taq-Probe polymerase<sup>c</sup></b>	Singleplex: 2U
	Multiplex up to four templates: 4U
RNA template <sup>d</sup>	As low as single digit copies of target RNA
Nuclease-free H <sub>2</sub> O	To a total volume of 20µl

**Footnotes of Table 4**

<sup>a</sup> The primer T<sub>m</sub> should be designed ≥60°C, preferably between 62°C to 65°C, using primer3 software for high efficiency and specificity.

<sup>b</sup> The probe's T<sub>m</sub> should be designed between 70-75°C.

<sup>c</sup> Taq-Probe polymerase, a specially engineered enzyme for S-shaped curve, is not included (Taq-Probe Polymerase, Cat W145-TP).

<sup>d</sup> RNA templates should be extracted by a qualified silica-based kit and eluted with low EDTA TE buffer (10mM Tris-HCl, 0.1mM EDTA, pH 8.0-8.3).

**Table 5. Compatible instruments**

RT-PCR Instrument	ROX required by instrument	Passive dye setup
Bio-Rad® iQ™5, CFX96, CFX384, Opticon Roche Lightcycler® Qiagen Rotor-Gene™ Eppendorf Mastecycler® Cepheid® SmartCycler®	Not recommended	Not necessary
Applied Biosystems® 7500, 7500 Fast, QuantStudio™, ViiA7™, Agilent Mx™	Low ROX (50nM final concentration)	Turn on ROX passive reference dye button
Applied Biosystems® 5700, 7000, 7300, 7700, 7900, 7900HT, 7900HT Fast, StepOne™, StepOnePlus™	High ROX (500nM final concentration)	Turn on ROX passive reference dye button

**Table 6. Standard thermocycling protocol**

Stage	Temperature	Period	Number of cycles
I	60°C	10min	1
II	95°C	2min	1
III	95°C	10sec	35-40
	60°C, signal acquisition	60sec	

**Table 7. Three-Step Thermocycling Protocol**

Stage	Temperature	Period	Number of cycles
I	60°C	10min	1
II	95°C	2min	1
III	95°C	10sec	35-40
	60°C	30sec	
	68-72°C, signal acquisition	30sec	

**Footnotes of Tables 6 and 7**

The three-step thermocycling protocol in Table 7 increases overall DNA polymerase activity by 50%, a more effective protocol than Table 6.

The primer concentration used in Tables 6 and 7 is typically 0.15-0.2uM.

**Related Products**

- Thermophilic Reverse Transcriptase, Cat. # W140
- Taq-Probe Polymerase, Cat.# W145
- Taq-Fast Polymerase, Cat.# W148
- 1-Step 2X Fast RT-PCR Master Mix-SYBR Green, Cat.# W147
- 1-Step 2X RT-PCR Master Mix-TaqMan Probe, Cat.# W143
- 1-Step 2X Multiplex RT-PCR Master Mix-TaqMan Probe, Cat. # W146
- 1-Step 2X Super Multiplex RT-PCR Master Mix-TaqMan Probe, Cat. # W149