

## Taq-Probe Polymerase

### Type and Catalog Number

Taq-Probe Polymerase

Cat. # W145, 250 reactions x 20 $\mu$ L (No ROX)

### Intended Use

- This DNA-dependent DNA polymerase is specially engineered for real-time PCR and RT-PCR amplification with TaqMan probe, singleplex or multiplex.
- It is used together with Thermophilic Reverse Transcriptase (Cat 140) for RT-PCR.

### Characteristics

- Taq-Probe polymerase is specially engineered for TaqMan probe, its increased 5'-3' exonuclease activity generating S-shaped curve.

**Table 1. Taq-Probe polymerase properties**

5'-3' polymerase activity	Yes
5'-3' exonuclease activity	Yes
3'-5' exonuclease activity	No
Reverse transcriptase activity	Neglectable
Incorporation of modified nucleotides	Yes, such as dUTP, fluorescence dye-labeled dNTPs
Terminal transferase activity	Minimal

**Table 2. PCR and RT-PCR optimal parameters**

Optimal RTase amount	Singleplex: 0.5-1U/20 $\mu$ l reaction
	Multiplex: 1-2U/20 $\mu$ l reaction
Optimal Taq-Probe polymerase amount	Singleplex: 2U/20 $\mu$ l reaction
	Multiplex up to four templates: 4U/20 $\mu$ l reaction
Optimal temperature	72-75 $^{\circ}$ C
Heat inactivation	>96 $^{\circ}$ C
dNTP concentration	Each 200 $\mu$ M
MgCl <sub>2</sub> concentration	3mM
Primer concentration	Each $\geq$ 0.15 $\mu$ M, depending on primer design and

	thermocycling
TaqMan probe	Each 0.15-0.25 $\mu$ M
Template	RNA: as low as single digit copies of target RNA
	Human genomic DNA: $\leq$ 1000ng/20 $\mu$ l reaction
Product size	Preferably 75-150bp

### Unit Definition

One unit of the enzyme catalyzes incorporation of 10nmol of deoxyribonucleotides into polynucleotide in 30 min.

### Production Source

E. coli strain

### Transportation and Storage

The kit can be stored at  $\leq$ -20 $^{\circ}$ C for 24 months. The product can be transported below 4  $^{\circ}$ C for up to 3 days.

**Table 3. Kit Contents**

Content	Amount
Taq-Probe polymerase, 25U/ $\mu$ l	500U/2500U
5x Taq-Probe buffer C	2x1mL/10mL
User manual	1

### Setup Reaction and Thermocycling

- Thaw 5x Buffer and other reaction components at room temperature, mix each component, centrifuge, and then place on ice.
- Determine the total volume for the number of reactions, add 5-10% extra volume, and prepare assay mix of all components except RNA template. Mix the assay mix, centrifuge, and then place on ice.
- Aliquot the assay mix into PCR tubes or plates.
- Add RNA template to PCR tubes or plate.
- Seal tubes with flat, optically transparent caps or seal plates with optically transparent film.
- Mix and then briefly centrifuge the tubes or plate.
- Program PCR instrument with indicated thermo-cycling protocol.
- Load PCR tubes or plates and start to run.
- Perform data analysis according to the PCR instrument instructions.

**Table 4. Set up a 20µl of reaction**

Content	Amount and final concentration
5x Buffer C	4µl
dNTPs	Each 200µM
Primers <sup>a</sup>	Each ≥0.15µM
TaqMan probe <sup>b</sup>	Each 0.15-0.25µM
Thermophilic Reverse Transcriptase <sup>c</sup>	Singleplex: 0.5-1U
	Multiplex: 1-2U
Taq-Probe polymerase	Singleplex: 2U
	Multiplex: 4U
Template	RNA: as low as single digit copies of target RNA
	Human genomic DNA: ≤1000ng/20µl reaction
H <sub>2</sub> O	To 20µl

**Footnotes of Table 4**

<sup>a</sup> The primer T<sub>m</sub> should be designed ≥60°C, preferably between 62°C to 65°C, using primer3 software for high efficiency and specificity.

<sup>b</sup> The probe's T<sub>m</sub> should be designed between 70-75°C.

<sup>c</sup> Thermophilic Reverse Transcriptase, Cat RT-PCR 140, is not included.

**Table 5. Compatible instruments**

RT-PCR Instrument	ROX required by instrument	Passive dye setup
Bio-Rad <sup>®</sup> iQ <sup>™</sup> 5, CFX96, CFX384, Opticon Roche Lightcycler <sup>®</sup> Qiagen Rotor-Gene <sup>™</sup> Eppendorf Mastercycler <sup>®</sup> Cepheid <sup>®</sup> SmartCycler <sup>®</sup>	Not recommended	Not necessary
Applied Biosystems <sup>®</sup> 7500, 7500 Fast, QuantStudio <sup>™</sup> , ViiA7 <sup>™</sup> , Agilent Mx <sup>™</sup>	Low ROX (50nM final concentration)	Turn on ROX passive reference dye button
Applied Biosystems <sup>®</sup> 5700, 7000, 7300, 7700,	High ROX (500nM final concentration)	Turn on ROX passive

7900, 7900HT, 7900HT Fast, StepOne <sup>™</sup> , StepOnePlus <sup>™</sup>		reference dye button
--	--	----------------------

**Table 6. Standard thermocycling protocol**

Stage	Temperature	Period	Number of cycles
I	60°C	10min	1
II	95°C	2min	1
III	95°C	10sec	35-40
	60°C, signal acquisition	60sec	

**Table 7. Three-Step Thermocycling Protocol**

Stage	Temperature	Period	Number of cycles
I	60°C	10min	1
II	95°C	2min	1
III	95°C	10sec	35-40
	60°C	30sec	
	68-72°C, signal acquisition	30sec	

**Footnotes of Tables 6 and 7**

The three-step thermocycling protocol in Table 7 increases overall DNA polymerase activity by 50%, a more effective protocol than Table 6.

The primer concentration used in Tables 6 and 7 is typically 0.15-0.2µM.

**Related Products**

- Thermophilic Reverse Transcriptase, Cat. # W140
- Taq-Probe Polymerase, Cat. # W145
- Taq-Fast Polymerase, Cat. # W148
- 1-Step 2X Fast RT-PCR Master Mix-SYBR Green, Cat. # W147
- 1-Step 2X RT-PCR Master Mix-TaqMan Probe, Cat. # W143
- 1-Step 2X Multiplex RT-PCR Master Mix-TaqMan Probe, Cat. # W146
- 1-Step 2X Super Multiplex RT-PCR Master Mix-TaqMan Probe, Cat. # W149