

## 1-Step 2X Fast RT-PCR Master Mix-SYBR Green Dye

### Product Name and Catalog Number

1-Step 2X Fast RT-PCR Master Mix-SYBR Green, Cat. W147-NR, W147-LR, W147-HR, 2x1mL for 200 reactions x 20 $\mu$ L (No ROX, Low ROX or High ROX)

### Intended Use

- The 1-Step Fast 2X RT-PCR Master Mix is used for real-time qualitative and quantitative RT-PCR amplifications with SYBR Green dye particularly for fast thermocycling.
- The master mix is a premixed, 2X concentrated solution that has all the components except for gene-specific primers and RNA template.

### Kit Characteristics

- The kit is designed for singleplex RT-PCR with SYBR Green dye.
- For the reverse transcription step, this kit uses a highly efficient Thermophilic Reverse Transcriptase (US patent pending), which is a thermophilic type A polymerase, with optimal temperatures of 60-62°C.
- The RTase is easily heat-inactivated at  $\geq 90^\circ\text{C}$  for 1min.
- The RTase efficiently synthesizes a complementary DNA strand on RNA template from a gene-specific primer,  $\leq 1$  unit per 20 $\mu$ L of reaction.
- The RTase reversely-transcribes single digit copies of target RNA molecules consistently.
- The kit also contains *Taq*-Fast DNA polymerase which extends more than 300 bases with short cycling program.
- The concentrations of the primers are variable depending on assay designs and thermocycling protocols (Table 1).
- The preferred PCR product size is  $\leq 150$ bp.
- The kit has three formulations of ROX, Low ROX or High ROX concentrations for your choice.

### Kit Contents

2X Master Mix (2x1mL for 200 reactions x 20 $\mu$ L)

### Transportation and Storage

The kit can be transported at below 4°C for up to 3 days.

The kit should be kept stable in the dark at -20°C for  $\leq 24$  months with  $\leq 10$  times of freeze-thaw cycles. The kit can be stored at 4°C for a week.

### Setup Reaction and Thermocycling

- Thaw 1-Step 2X RT-PCR Master Mix and other reaction components at room temperature, mix each component, centrifuge, and then place on ice.
- Determine the total volume for the number of reactions, add 5-10% extra volume, and prepare assay mix of all components except RNA template. Mix the assay mix, centrifuge, and then place on ice.
- Aliquot the assay mix into PCR tubes or plates.
- Add RNA template to PCR tubes or plate.
- Seal tubes with flat, optically transparent caps or seal plates with optically transparent film.
- Mix and then briefly centrifuge the tubes or plate.
- Program PCR instrument with indicated thermo-cycling protocol.
- Load PCR tubes or plates and start to run.
- Perform data analysis according to the PCR instrument instructions.

**Table 1. Setting up a 20 $\mu$ L or 10 $\mu$ L reaction**

Component	Volume per 20 $\mu$ L	Volume per 10 $\mu$ L	Final concentration
2X Master Mix	10 $\mu$ L	5 $\mu$ L	1X
Primers <sup>a</sup>	Variable	Variable	Each 150-900nM
RNA template <sup>b</sup>	Variable	Variable	As low as single digit copies of target RNA to $\leq 1\mu\text{g}$ total RNA
H <sub>2</sub> O	To 20 $\mu$ L	To 10 $\mu$ L	

### Footnotes of Table 1

<sup>a</sup>The primer's  $T_m$  should be designed  $\geq 60^\circ\text{C}$ , preferably between 62°C to 65°C, using primer3 software for high efficiency and specificity.

<sup>b</sup>RNA templates should be extracted by a qualified silica-based kit and eluted with low EDTA TE buffer (10mM Tris-HCl, 0.1mM EDTA, pH 8.0-8.3).

**Table 2. Compatible instruments**

RT-PCR Instrument	ROX required by instrument	Passive dye setup
Bio-Rad <sup>®</sup> iQ <sup>™</sup> 5, CFX96, CFX384, Opticon Roche Lightcycler <sup>®</sup> Qiagen Rotor-Gene <sup>™</sup>	Not recommended	Not necessary

Eppendorf Mastercycler® Cepheid® SmartCycler®		
Applied Biosystems® 7500, 7500 Fast, QuantStudio™, ViiA7™, Agilent Mx™	Low ROX (50nM final concentration)	Turn on ROX passive reference dye button
Applied Biosystems® 5700, 7000, 7300, 7700, 7900, 7900HT, 7900HT Fast, StepOne™, StepOnePlus™	High ROX (50nM final concentration)	Turn on ROX passive reference dye button

**Table 3. Standard thermocycling protocol**

Stage	Temperature	Period	Number of cycles
I	60°C	10min	1
II	95°C	2min	1
III	95°C	12sec	35-40
	60°C, signal acquisition	60sec	
IV	60°C to 95°C	Various	1

**Table 4. Three-Step Thermocycling Protocol**

Stage	Temperature	Period	Number of cycles
I	60°C	10min	1
II	95°C	2min	1
III	95°C	10sec	35-40
	60°C	30sec	
	68-72°C, signal acquisition	30sec	
IV	68°C to 95°C	Various	1

**Footnotes of Tables 3 and 4**

The three-step thermocycling protocol in Table 4 increases overall DNA polymerase activity by 50%, a more effective protocol than Table 3.

The primer concentration used in Tables 3 and 4 is typically 0.15-0.2uM.

**Table 5. Fast thermocycling protocol**

Stage	Temperature	Period	Number of cycles
I	60°C	5-10min	1
II	95°C	1min	1
III	95°C	5sec	35-40
	60°C, signal acquisition	30sec	
IV	60°C to 95°C	Various	1

**Footnotes of Table 5**

The product size for the fast-thermocycling protocol is preferred to be less than 90bp.

The primer concentration used is typically between 0.4uM and 0.9uM.

**Related Products**

- Thermophilic Reverse Transcriptase, Cat. W140
- Taq-Probe Polymerase, Cat. W145
- Taq-Fast Polymerase, Cat. W148
- 1-Step 2X Fast RT-PCR Master Mix-SYBR Green, Cat. W147
- 1-Step 2X RT-PCR Master Mix-TaqMan probe, Cat. W143
- 1-Step Multiplex 2X RT-PCR Master Mix-TaqMan probe, Cat. W146
- 1-Step 2X Super Multiplex RT-PCR Master Mix-TaqMan Probe, Cat. W149

**Precautions**

If you order a “No ROX” master mix but have an Applied Biosystems/ThermoFisher instrument, please **turn off ROX passive reference dye button** when setup assays.