

## 2X Fast qPCR Master Mix-SYBR Green

### Product Name and Catalog Number

2X Fast qPCR Master Mix-SYBR Green  
Cat. # W157-NR, W157-LR, W157-HR,  
Size: 2 x 1ml for 200 reactions (No ROX, Low  
ROX or High ROX)

### Intended Use

- This 2X Fast qPCR Master Mix is used for real-time qualitative and quantitative qPCR with SYBR Green dye. This product is for research use only.
- The master mix is a premixed, 2X concentrated solution that has all the components except for gene-specific primers and DNA template.

### Kit Features

- The kit is designed for singleplex qPCR with SYBR Green dye.
- The kit contains *Taq*-Fast DNA polymerase which can extend more than 300 bases with a short cycling program.
- The concentrations of the primers are variable depending on assay designs and thermo-cycling protocols (Table 1).
- The preferable PCR product size is  $\leq 150$ bp.
- The kit has three formulations of ROX, Low ROX or High ROX concentrations for your choice (see Table 2).

### Kit Contents

2X Master Mix (2x1ml for 200 reactions)

### Transportation and Storage

The kit can be transported at  $\leq 4^{\circ}\text{C}$  for up to 3 days.

The kit should be stored in the dark at  $-20^{\circ}\text{C}$  for no more than one year with no more than 10 times of freeze-thaw cycles. The kit can be stored at  $\leq 4^{\circ}\text{C}$  for a week.

**Table 1. Setting Up a 20 $\mu\text{L}$  or 10 $\mu\text{L}$  reaction**

Component	Volume per 20 $\mu\text{L}$	Volume per 10 $\mu\text{L}$	Final concentration
2X Master Mix	10 $\mu\text{L}$	5 $\mu\text{L}$	1X
Primers <sup>a</sup>	Variable	Variable	Each 150-900nM
DNA template <sup>b</sup>	Variable	Variable	$\leq 60$ ng or $\leq 30$ ng genomic DNA
H <sub>2</sub> O	To 20 $\mu\text{L}$	To 10 $\mu\text{L}$	

### Footnotes of Table 1

<sup>a</sup> Each primer's  $T_m$  should be designed  $\geq 60^{\circ}\text{C}$ , preferably between  $62^{\circ}\text{C}$  to  $65^{\circ}\text{C}$ , using primer3 software for high efficiency and specificity.

<sup>b</sup> DNA templates should be extracted by a qualified silica-based kit and eluted with low EDTA TE buffer (10mM Tris-HCl, 0.1mM EDTA, pH 8.0-8.3).

### Applicable Instruments

**Table 2. Compatible instruments**

PCR Instrument	ROX required by instrument	Passive dye setup
Bio-Rad <sup>®</sup> iQ <sup>™</sup> 5, CFX96, CFX384, Opticon Roche Lightcycler <sup>®</sup> Qiagen Rotor-Gene <sup>™</sup> Eppendorf Mastercycler <sup>®</sup> Cepheid <sup>®</sup> SmartCycler <sup>®</sup>	Not recommended	Not necessary
Applied Biosystems <sup>®</sup> 7500, 7500 Fast, QuantStudio <sup>™</sup> , ViiA7 <sup>™</sup> , Agilent Mx <sup>™</sup>	Low ROX (50nM final concentration)	Turn on ROX passive reference dye button
Applied Biosystems <sup>®</sup> 5700, 7000, 7300, 7700, 7900, 7900HT, 7900HT Fast, StepOne <sup>™</sup> , StepOnePlus <sup>™</sup>	High ROX (500nM final concentration)	Turn on ROX passive reference dye button

## Setting Up Thermal Cycling

**Table 3. Standard thermo-cycling Protocol**

Stage	Temperature	Period	Number of cycles
I	95°C	2min	1
II	95°C	12sec	35-40
	60°C, signal acquisition	60sec	
III	60°C to 95°C	Various	1

### Footnotes of Table 3

The primer concentration used is typically 0.2uM.

**Table 4. Fast thermocycling protocol**

Stage	Temperature	Period	Number of cycles
I	95°C	1min	1
II	95°C	5sec	35-40
	60°C, signal acquisition	30sec	
III	60°C to 95°C	Various	1

### Footnotes of Table 4

The product size for the fast thermocycling protocol is preferred to be less than 90bp.

The primer concentration used is typically between 0.4uM and 0.9uM.

## Quality Control

Not detectable DNase and RNase contaminations.

## Related Products

- 2X qPCR Master Mix-TaqMan probe, Cat.# W153 (No ROX, High ROX or Low ROX)
- 2X Multiplex qPCR Master Mix-TaqMan probe, Cat.# W156 (No ROX, High ROX or Low ROX)

## Precautions

If you order a “No ROX” master mix but you have an Applied Biosystems/ThermoFisher instrument, please **turn off ROX passive reference dye button** when setup assays.