

## 101Transfection Kit

**Cat. #:** [P901](#)

**Storage:** 4°C.

**Shelf life:** 6 months

**Product Size:** 1 mL 101Transfection Reagent

**Product Description** (This product is for research use only)

101Transfection Reagent is a **unique nanoparticles** that delivers more DNA and siRNA to cells than other lipid-based transfection kits.

- ✓ **No media changes** needed
- ✓ **Antibiotics and serum compatible** – no need to remove
- ✓ **One-step** incubation - 15 minutes
- ✓ Work on **most cell types**

**Reagent required not provided in this kit:** serum-free DMEM with High Glucose

**Positive control:** human embryonic kidney 293T cells

**Table 1. Suggested Reaction Matrix**

Culture Dish	Area (cm <sup>2</sup> )	Cell Number	Medium (mL)	DNA (µg)	101Transfection Reagent (µL)	DMEM (µL)
96-Well	0.2	1.4 x 10 <sup>4</sup>	0.1	0.1	0.3	10
48-Well	1	9.5 x 10 <sup>4</sup>	0.25	0.25	0.75	20
24-Well	2	1.9 x 10 <sup>5</sup>	0.5	0.5	1.5	50
12-Well	4	3.8 x 10 <sup>5</sup>	1.0	1	3	100
6-Well /35 mm	10	9.5 x 10 <sup>5</sup>	2.0	2.5	7.55	200
60 mm / T25	25	2.1 x 10 <sup>6</sup>	5.0	6-8	15-24	300
100 mm / T75	75	5.5 x 10 <sup>6</sup>	10	15-20	35-45	500
150 mm /T150	150	15 x 10 <sup>6</sup>	20	25-30	60-80	1000

**Note:**

For different cell type, the seeding cell density may vary. Customer should test the seeding condition to ensure the 70 ~ 80% confluency at the time of transfection. The transfection efficiency also depend on the cell type. For difficult-to-be-transfected cells, customer can try to increase the 101Transfection reagent volume up to 4 µl per µg of DNA.

**Protocol** (example of transfecting 293T cells in **one well of 24-well plate**. Refer to Table 1 for other plates or dishes)

**Notice:** Transfect cells at **70-80% confluency** to achieve high transfection efficiency and low toxicity.

1. **Plating:** seed the cells **18 to 24 hours** prior to transfection in 24-well plate, at the density that allow the cells reaches 70~80% confluency at the time of transfection.
2. **Change medium: 2 hours before** transfection, remove culture medium and add **0.5 mL fresh** complete culture medium. (No need to use serum-free / antibiotics-free medium.)
3. In **tube 1** add: **0.5 µg DNA**  
**50 µl DMEM** (serum-free, High Glucose)  
Pipet up and down to mix well
4. In **tube 2** Add: **1.5 µl 101Transfection Reagent**  
**50 µl DMEM** (serum-free, High Glucose)  
Gently pipet up and down to mix well
5. Add **tube 2** into **tube 1**. **Vortex for 10 seconds**.
6. Incubate for **15 minutes** at **room temperature**.
7. Add the incubated mixture **drop-wise** to the cells, and homogenize by **gently swirling** the plate.
8. Return the cells to incubator.
9. Check transfection efficiency 24 to 48 hours post transfection.

-- The end --