

## Thermophilic Reverse Transcriptase

**Product Name:** Thermophilic Reverse Transcriptase:

**Cat. #:** W140 (200 unit / reaction)

### Applications

- The reverse transcriptase is a thermophilic type A polymerase (US patent pending) with optimal temperatures at 60-62°C, and can be heat-inactivated at ≥90°C.
- The enzyme efficiently synthesizes a complementary DNA strand initiating from a gene-specific primer, *one unit per 20 µL of reaction*.
- The enzyme can detect single digit copies of target RNA depending on assay design and optimization.
- It is particularly suitable to one-step real-time qualitative and quantitative RT-PCR.

### Enzyme properties

Reverse transcriptase activity	Yes
5'-3' DNA polymerase activity	Yes
5'-3' exonuclease activity	Yes
3'-5' exonuclease activity	No
Incorporation of modified nucleotides	Yes, such as dUTP, fluorescence dye-labeled dNTPs
Heat-Inactivated	≥90°C for 2 min
DNase or RNase activity	Not detectable

### Optimal parameters

Optimal RTase amount	1-2 U per 20 ul reaction
Optimal temperature	60-62°C
Optimal incubation time	10 min
Mg <sub>2</sub> concentration	1.5-3 mM
dNTP concentration	Each 40-200 uM
Product size	Preferred ≤150 bp

### Unit definition

One unit of the enzyme catalyzes the incorporation of 1 nmol of deoxyribonucleotides into polynucleotide in 10 min.

### Production source

*E.coli* strain

Product Component	Amount
RTase, 10 U/µl	200 Units
10x RT-PCR buffer-SYBR Green dye: 500mM Tris-HCl (pH 8.3 @ 25°C), 150mM (NH <sub>4</sub> ) <sub>2</sub> SO <sub>4</sub> , 15mM MgCl <sub>2</sub> , 0.2% Tween20, 1mM DTT, 500 µg/ml BSA, with SYBR Green dye	1,000 µl
10x RT-PCR buffer-SYBR Green dye: 500mM Tris-HCl (pH 8.3 @ 25°C), 150mM (NH <sub>4</sub> ) <sub>2</sub> SO <sub>4</sub> , 15mM MgCl <sub>2</sub> , 0.2% Tween20, 1mM DTT, 500 µg/ml BSA, with SYBR Green dye	1,000 µl

**Notes:** **No ROX passive reference dye** is in the 10x RT-PCR buffers.

### Ship and Store

The kit can be shipped at **4°C** (for up to 3 days).

It should be stored at **-20°C** for 24 months before use.

### Set up 20µl of reaction

Component	Amount or final concentration
10x RT-PCR buffer	2 ul
dNTPs	Up to 200 uM
Target specific primers <sup>a</sup>	Each ≥0.1 uM
Taqman probe or SYBR Green <sup>b</sup>	Variable
RTase	1-2 U
<i>Taq</i> DNA polymerase <sup>c</sup>	1-2 U
RNA template <sup>d</sup>	As low as single digit copies of target RNA
Nuclease-free H <sub>2</sub> O	To a total volume of 20 µl

<sup>a</sup> The primer T<sub>m</sub> should be designed ≥60°C using primer3 software for high efficiency at the optimal temperature.

<sup>b</sup> TaqMan probe or SYBR Green dye can be used for fluorescent signal.

<sup>c</sup> Not included in the kit.

<sup>d</sup> RNA templates should be extracted by a qualified silica-based kit and eluted with low EDTA TE buffer (10 mM Tris-HCl, 0.1 mM EDTA, pH 8.0-8.3).

### Suggested thermo-cycling for RTase

Start reverse transcription at 60°C for 10 minutes, and then inactivate Thermophilic Reverse Transcriptase at 94-95°C for 2 minutes, followed by a PCR program.

**Note: Turn off ROX passive reference dye button** when setup assays on Applied Biosystems/ThermoFisher instruments.