

## Magnetic Beads (DNA)

**Cat. #:** P920-1 (1 mL); P920-5 (5 mL); P920-30 (30 mL); P920-60 (60 mL); P920-450 (450 mL)

**Storage:** 4°C (Do not freeze, keep sterile)

**Shelf Life:** 12 months

### Product Description:

For DNAs (dsDNA or ssDNA,  $\geq 100$ base) purification in NGS library preparation, microarray sample preparation, PCR products purification, Sanger sequencing, gene cloning, and many other procedures with DNA purifications. This product is for research use only.

\* Additional magnetic racks for tubes or plates are required.

## Protocol

**Important note:** Precipitates during storage is normal. **Before each use**, mix well by shaking, thoroughly pipetting or vortexing to homogenous solution. No precipitate should be seen at the tube bottom. Well mixing is critical for high yield of DNA purification.

### Magnetic Beads Prepare

1. Before each use, shake the beads bottle well to homogenous solution. Transfer needed volume of the Magnetic Beads solution to a fresh tube, and **warm it up** at room temperature for 15 minutes.

### Add Magnetic Beads to DNA sample

2. Vortex thoroughly again to mix well. In one 1.5ml Eppendorf tube or 96-well plate, add **1.8 volume** of beads solution to 1 volume of DNA sample solution, mix the sample/bead mixture by gently pipetting for 15 times or vortexing for 2 minutes.

**Note:** for best result, prepare DNA concentration at **10 ng/uL ~ 500 ng/uL**, and in **10 uL ~ 100 uL** volume.

### Incubate at room temperature

3. Let the mixture stand at **room temperature** for **5 minutes**. Briefly spin-down the mixture in microcentrifuge, and load it to a **magnetic rack** (e.g., Life Technologies MagnaRack Cat. No. CS15000 or Magnetic Stand-96 Cat. No. AM10027).

### Removal of the Buffer

4. Let the mixture stand at **room temperature** for **3 minutes**. The solution will become clear after the **integrated sample/bead unites** are collected at the tube wall.
5. Remove the liquid phase from the tube without disturbing the integrated sample/bead unites on the tube wall.

### Wash

6. Keep the tube on the magnetic rack, add **500uL fresh 70% ethanol/isopropanol** to the pellet and wash for 30 seconds (If the wash within a tube, turn the tube around for 180° on the magnetic rack for better wash.).
7. Remove the liquid with a pipette.
8. Repeat steps 6 and 7 once.
9. Spin the tubes or plates at **12,000 rpm** for **20 seconds** and reload to the magnetic rack for **10 seconds**.
10. Use a 10uL pipetting tips to remove the leftover liquid inside the tube as complete as possible, without disturbing the integrated sample/bead unites

### Air Dry

11. Let the tube cap / plate cap open for **5 minutes** to air dry.

### Elution

12. Add 20 ul ~ 40 uL TE or ddH<sub>2</sub>O to the pellet and pipet to mix well. Incubate at room temperature for 2 minutes.
13. **Briefly spin** the tube / plate, and reload it back to the **magnetic rack** for **2 minutes**.
14. Transfer the **supernatant (eluted DNA)** to a new Eppendorf tube or plate.
15. The purified sample is ready for use, or store at 4°C for short term or -20°C for long term.

-- The end --