

All-Lab Water Bath Treatment Tablet (1 tablet / treatment of 1 – 2 Gal.)

Cat. #: T20 (20 Tablets), T200 (200 Tablets)

Shipping and Storage: Ship and stored at room temperature. Store in cool, dry and dark place.

Shelf Life: 12 months

Product Description:

This product can effectively kill all bacteria and fungus commonly found in water bath and incubator water tray. Reduces hard water deposits to extend water bath life span. Easy to use in tablet form (1 tablet treats 1-2 gallons of water, max. 4 gallons). After the tablet is dissolved in 1-2 gallons of water, the water is Not harmful to skin. Twice a week is recommended to keep your water bath clean.

Product Contents

Component	T20	T200
All-Lab Water Bath Treatment Tablet	20 tablets	200 tablets

Direction of Use:

- 1) Put one tablet in the water bath or incubator water tray of 1 – 2 gallons water, max. 4 gallons, circulate the water for 1 minute to dissolve it evenly in water.
- 2) After 10 minutes, drain the whole water, and remove all debris if any from the water container. Thoroughly clean the container with clean water, done. One tablet treatment is good for one week. This product is for research use only. Follow all related State requirements concerning cleaning if any.

STORAGE AND DISPOSAL

Do not contaminate drink water, food, or feed during storage or disposal. Store in original container in areas inaccessible to children. Store at room temperature lower than 100 F. Do not store in direct sunlight. Do not store on side. Avoid creasing or impacting of side walls. Discard in trash or recycle. Wastes resulting from the use of this product may be disposed of on site or at an approved waste disposal facility.

Active Ingredients

1-bromo-3-chloro-5,5-dimethylhydantoin	3.6%	1,3-dichloro-5,5-dimethylhydantoin	1.9%
1,3-dichloro-5-ethyl-5-methylhydantoin	1.1%	Sodium bromide	4.9%
Inert Ingredients:	<u>88.5%</u>		
Total:100.0%			

PRECAUTIONARY STATEMENTS : HAZARDS TO HUMANS AND DOMESTIC ANIMALS. KEEP OUT OF REACH OF CHILDREN. After the tablet is dissolved in water, the water is not harmful to skin.

WARNING: Before the tablet is dissolved in water, it causes substantial but temporary eye injury. Harmful if

absorbed through skin, eyes or clothing. Do not get in eyes, on skin or on clothing. Wear goggles, face shield or safety glasses. Wash thoroughly with soap and water after handling and before eating, drinking, chewing gum, using tobacco, or using the toilet. Remove and wash contaminated clothing before reuse.

FIRST AID: Take this manual with you when calling a Poison Control Center or doctor or going for treatment.

IF IN EYES: Hold eye open and rinse slowly and gently with water for 15 – 20 minutes. Remove contact lenses if present, after the first 5 minutes, then continue rinsing eye. Call a poison control center or doctor for treatment advice.

IF ON SKIN OR CLOTHING: Take off contaminated clothing. Rinse skin immediately with plenty of water for 15-20 minutes. Call a poison control center or doctor for treatment advice.

PHYSICAL AND CHEMICAL HAZARDS: STRONG OXIDIZING AGENT. Do not mix this product in concentrated form with any other chemicals. Mix only with water. Use clean dry utensils. Do not add this product to any dispensing device containing remnants of any other product. Such use may cause a violent reaction leading to fire or explosion. Contamination with moisture, organic matter, or other chemicals may start a chemical reaction with generation of heat, liberation of hazardous gases, and possible generation of fire and explosion. In case of contamination or decomposition, do not reseal container. If possible, isolate container in open air or well-ventilated area. Flood with large volumes of water, if necessary.
EPA Reg. No. 75562-1 EPA Est. No. 85642-NJ-001

Customers are also interested in

<i>Better than Matrigel</i>	Real 3D cell gel, 3-stiffness for different type cells # P720
<i>Mouse Tail DNA Extract</i>	20 min. fast. Accurate & consistent. No use of organic solvents. Ready for PCR Cat.# T605
<i>1 drop do PCR Mix</i>	PCR Master Mix - easy dripping Cat.# W2599-5
<i>10x Virus Titer-up Boost</i>	Increases Lenti/Retro virus titer up to 10-fold # P906/ P909
<i>Lentiviral Expression Vectors</i>	4 promoters: SFFV, CAG, CMV, EF1. 4 selections
<i>2x Flash HS PCR Mix (Dye)</i>	Fast extension rate (5-15 sec/kb) and fast HotStart # W3007-1
<i>Simple and Reliable</i>	SV40T / hTERT <i>Cell Immortalization Kits</i>