

Ultra SYBR Green qPCR Master Mix (2X, with ROX I)

Cat. #: W2601-1 (1 mL); W2601-5 (5 mL)

Storage: Store at -20°C for up to 1 year and avoid freeze-thaw cycles. Stored at 4°C if it is frequent used within short time.

Platform: ABI Prism7500/7500 Fast, Stratagene Mx3000/Mx3005P, Corbett Rotor Gene 3000

Product Components

Cat. No.	Kit Size	2x Ultra SYBR Green qPCR Master Mix (with ROX I)	RNase-Free Water
W2601-1	50 µL volume x 40 reactions,	1 mL	1 mL
W2601-5	50 µL volume x 200 reactions	1 mL x 5	5 mL

Application

- Gene expression analysis
- Microarray result validation

Product Description (This product is for research use only.)

UltraSYBR Mixture (With ROX) provides a convenient real-time PCR specific reagent by intercalater method using SYBR Green I for detection. It is a 2x concentrated premix of GoldStar Taq DNA polymerase, dNTPs, Mg²⁺, SYBR Green I dye (detection), ROX reference dye and PCR buffer components. The product utilizes an enzyme for hot start, GoldStar Taq DNA polymerase, which confers a significant reduction in non-specific PCR amplification. As the enzyme buffer system is optimized for real time PCR, the product offers high amplification efficiency and high detection sensitivity in real-time PCR.

In light of the fact that the qPCR instruments can vary from user to user, Cowin offers the UltraSYBR Mixture in a range of formulations, each of which has been carefully optimized to confer the best performance according to the make and model of a qPCR machine. Please use the following table as a guide for selecting the UltraSYBR Mixture that will be most compatible with your choice of a particular instrument/model.

Key features

- Easy-to-use 2X premix reagent including SYBR Green I
- Hot start PCR enzyme enables high efficiency and sensitivity in real-time PCR (qPCR)
- Accurate detection and quantification of target gene through real time PCR.

Precautions before Use

- Prior to use, make sure the reagent is evenly mixed by inverting the bottle several times without creating bubbles and briefly centrifuge the tube before use.
- Store at 2-8°C to avoid repeated freeze-thaw cycles if it is frequent used within short time.
- This product contains SYBR Green I and ROX Reference Dye. Avoid exposing to strong light when preparing the reaction mixture.
- Use fresh disposable tips to avoid contamination between samples when preparing or dispensing reaction mixtures.

Protocol

The following example is applied to conventional real-time PCR reaction system. The reaction conditions can be optimized according to the template, primer structure. This protocol is for a reaction size of approximately 50 μ L. The reaction size may be adjusted as desired.

1. Prepare the PCR mixture shown below

Reagent	50 μ L PCR reaction	Final Concentration
2x Ultra SYBR Green qPCR Master Mix (with ROX I)	25 μ L	1x
Forward Primer, 10 μ M	1 μ L	0.2 μ M
Reverse Primer, 10 μ M	1 μ L	0.2 μ M
Template DNA	2 μ L	-
RNase-Free Water	Up to 50 μ L	-

Note:

- 1) A final primer concentration of 0.2 μ M is likely to yield good results. However, if there is an issue with reactivity, use a primer concentration between 0.1 and 1.0 μ M.
 - 2) It is preferable to use 10-100 ng genomic DNA or 1-10 ng cDNA as template. The optimal quantity varies depending on the number of target copies in the template solution. Make serial dilutions to determine the appropriate amount.
 - 3) If cDNA (RT reaction mixture) is used as template, the template volume should be no more than 10% of the PCR mixture.
2. Briefly centrifuge reaction tubes then set them for real-time PCR reaction.

This protocol is set according to ABI 7500 Real-Time PCR instrument. Two-step PCR reaction procedure is recommended. Try this protocol first and optimize PCR conditions as necessary.

Please notice that the initial denaturation step prior to PCR should be at 95°C for 10 min)

Two-step PCR protocol

Procedure	Temperature	Time	
Initial denaturation	95°C	10 min	
Denaturation	95°C	15 s	} 35-40 cycles
Annealing / Extension	60°C	1 min	
Melting curve analysis			
	95°C	15 s	
	60°C	1 min	
	95°C	15 s	
	60°C	15 s	

Note:

- 1) This product utilizes Goldstar Taq polymerase, which is an enzyme for hot start PCR. Initial denaturation step prior to PCR should be at 95°C for 10 min
- 2) This protocol is set according to ABI 7500 Real-Time PCR, melting curve analysis protocol should be adjusted according to the qPCR instrument.
- 3) Please try annealing of two-step PCR at 60°C for 30 seconds at first. The temperature should be optimized within the range of 60 - 64°C if optimization is required.
- 4) Perform a three-step PCR when using primers with low T_m values. The annealing temperature for three-step PCR should be optimized within the range of 56-64°C.

-- The end --