

SARS-CoV-2 Pseudovirus Infection Enhancer (CoV-2-PIE)

Cat. #: CoV2-1ml, CoV2-200 works perfectly for Luc Pseudovirus to get robust signal.

This viral infection enhancer specifically designed to facilitate coronavirus and pseudovirus infection of target cells. It is designed to penetrate the cortical actin barrier, thereby 5 to 20 folds enhancing productive viral infection of a variety of host cells by coronaviruses and pseudoviruses. **CoV-2-PIE** is based on the scientific theory that the actin cytoskeleton is a natural barrier for viral entry and post-entry intracellular migration (Yoder et al., *Cell*, 2008, 134:782). It's formulated to mainly enhance viral infection of adherent target cells, which are frequently also limited by poor absorption of virion particles to target cells. CoV-2-PIE is intended for Research Use. [Link to web page](#)


Applications

- Enhancing coronavirus and pseudovirus transduction of target cells
- Enhancing infection of other enveloped viruses
- Facilitating recovery of infectious viruses from cell or tissue cultures
- Facilitating anti-viral drug and neutralizing antibody screening efficiency

Important notes

- 1) **CoV-2 PIE** is formulated as 10 X concentrated, **works perfectly for Luc pseudovirus**. Once thawed, **CoV-2-PIE** should be stored at 4°C, stable for 6 months. *Do not re-freeze or leave it at room temperature.*
- 2) This viral infection enhancer works with most cell lines to enhance viral infection. On average, **CoV-2-PIE** enhances productive viral infection by **5- to 20-folds**. (depends on *the types of viruses and cells.*)

**SARS-CoV-2
pseudovirus infection of Vero E6
(with infection enhancer, Cov2-1ml)**



Example of results from our customers:

SARS-CoV-2-S(Luc) pseudovirus infection of Vero E6 cells (with pseudovirus infection enhancer).

Vero E6 cells were transduced with SARS-CoV-2-S(Luc) lenti-pseudovirus (with a luciferase reporter) in the presence or absence of infection enhancer. Reporter expression was quantified at 3 days post transduction (luciferase assay).

Protocol

Example - CoV-2-PIE enhances SARS-CoV-2 lenti-pseudovirus transduction of Vero E6 cells:

- 1) Count Vero E6 cells to be infected and seed $\sim 1 \times 10^5$ cells per well into 12-well plates (0.5 ml per well). Culture cells until cells stably adhere to the plates (4–12 hours).

Note: Cell viability should be $\geq 80\%$.

- 2) Before infection, wash cells with 2 ml medium, and leave 250 μ l medium in each well.
- 3) Pre-treat cells by adding 25 μ l of **CoV-2-PIE (10 X)** so that the **CoV-2-PIE** concentration is **1 X**. Mix and incubate for 30–60 minutes.
- 4) Add virus to the cells and mix. **Note volume of virus used.**
- 5) Add **CoV-2-PIE (10 X)** in an amount equal to 1/10 of the virus volume used, *e.g.*, if 100 μ l of virus is used, add 10 μ l of **CoV-2-PIE**. Incubate at 37°C for 4-6 hours.
- 6) Add 2 ml fresh media to wash cells.
- 7) After washing, add 2 ml fresh complete medium.
- 8) Culture infected cells for 2-3 days to signal detection.

Catalog No.	Product Size
CoV2-200	200 μ l of 10 X CoV-2-PIE solution
CoV2-1ml	1 ml of 10 X CoV-2-PIE solution

We also provide

Ready to use SARS-CoV and SARS-CoV-2 S protein pseudotyped lentiviruses. [Link](#)

HIV Rev-dependent Reporter Cells (Cat.# HRC-1) and HIV infection enhancer (H901-1) [Link](#)

Real 3D Cell Culture Gel Col-Tgel (Cat.# P720) [Link](#)

Lentivirus / Retrovirus 10X Titer-Up Reagent (Cat.# P906 / P909) [Link](#)

95% pure Exosome Isolation kits [Link](#)

Mouse tail DNA extraction kit [Link](#)