

## 293T Transfection Kit for **Lentivirus Packaging**

**Cat. #:** [P902](#)

**Storage:** 4°C.

**Shelf life:** 6 months

**Product Size:** 1 mL 293T Transfection Reagent for Lentivirus Packaging

**Product Description** (This product is for research use only)

This transfection reagent is a **unique recipe** for delivery of Lentivirus packaging DNA to 293T cells at highest efficiency, to produce high titer lentivirus.

- ✓ Higher efficiency than lipid-based transfection Reagents
- ✓ **No media changes** needed
- ✓ **Serum compatible**
- ✓ **One-step** incubation - 15 minutes

**Reagent required not provided in this kit:** serum-free DMEM with High Glucose

**Table 1. Suggested Reaction Matrix**

Culture Dish	Area (cm <sup>2</sup> )	Cell Number	Medium (mL)	DNA (µg)	101Transfection Reagent (µL)	DMEM (µL)
6 - Well	10	9.5 x 10 <sup>5</sup>	2.5	4	7.55	200
60 mm	25	2.5 x 10 <sup>6</sup>	4.0	7-12	15-24	300
100 mm	55	5.5 x 10 <sup>6</sup>	10	15-25	35-45	500
150 mm	150	15 x 10 <sup>6</sup>	20	45-75	60-80	1000

**Protocol** (example of producing virus by transfecting 293T cells in **one 100mm petri dish**.)

Refer to Table 1 for other plates or dishes)

**Notice:** Transfect cells at **70-90% confluency** for high transfection efficiency and low toxicity. Before plating, pipet the 293T cell suspension vigorously to achieve **single cells**. Plate the cells in **evenly** distribution manor.

### Day 1: Plating cells

1. **Plating: 18 to 24 hours** prior to transfection, coat 100mm dishes with **6 mL 1x Gelatin** (101Bio, Cat. #: P910) for 60 min. Aspirate gelatin, and plate **5.0 x 10<sup>6</sup>** fast growing 293T cells.

## Day 2: Transfection

- Change medium: 2 hours before** transfection, remove culture medium and add **10 mL fresh** complete culture medium (with 10% serum / **without antibiotics**)
- In **tube 1** add:  
5 µg Lentivirus vector DNA (Containing your gene of interest)  
20 µg Lentivirus Packaging Mix (101Bio, Cat. # P904P)  
500 µL DMEM (serum-free, High Glucose)  
Pipet up and down to mix well
- In **tube 2** Add:  
40 µL 293T Transfection Reagent for Lentrivirus Packaging  
500 µL DMEM (serum-free, High Glucose)  
Gently pipet up and down to mix well
- Incubate at room temperature (**20–25°C**) for **3 min**.
- Add **tube 2** into **tube 1**, pipet up and down several times. **Vortex** for **10 seconds**.
- Incubate for **15 minutes** at **room temperature**.
- Add the incubated mixture **drop-wise** to the cells, and **gently swirl** the plate to disperse evenly in the plate.
- Return the cells to 37°C incubator with 5% CO<sub>2</sub>.

## Day 3: Add Titer-Up to increase virus titer

Add **20 µl** of Lentivirus **Titer-Up** (101Bio, Cat. #: 909) to the medium. Return the plates to the cell culture incubator.

## Day 4: Collect virus

- Collect virus supernatant twice at **48 and 72 hours** post transfection into a 50mL sterile FALCON tube. Centrifuge at **3,000rpm** for **15 minutes** at 4°C to remove cell debris. Filter the clear supernatant through **0.45 µm** syringe tip filter.
- The filtered clear supernatant is virus soup. Use it immediately or aliquot into sterile 1.5-mL tubes and store at **-80°C**, for **up to 3 months**.
- (Optional) To concentrate virus, add **5X Lentivirus Concentration Solution** (101Bio, Cat. #: P904C) to the viral supernatant (volume of Lentivirus Concentration Solution vs. volume of viral supernatant = 1:4) and mix thoroughly. Put the mixture to 4°C refrigerator overnight and spin the virus pellet down next day. Please refer our P904C user manual for details.

-- The end --