

293T Transfection Kit for **Retrovirus Packaging**

Cat. #: [P903](#)

Storage: 4°C

Shelf life: 6 months

Product Size: 1 mL 293T Transfection Reagent for Retrovirus Packaging

Product Description (This product is for research use only)

This transfection reagent is a **unique recipe** for delivery of **Retrovirus packaging DNA** to 293T cells at highest efficiency, to produce high titer retrovirus.

- ✓ Higher efficiency than lipid-based transfection Reagents
- ✓ **No media changes** needed
- ✓ **Serum compatible**
- ✓ **One-step** incubation - 15 minutes

Reagent required not provided in this kit: serum-free DMEM with High Glucose

Table 1. Suggested Reaction Matrix

| Culture Dish | Area (cm ²) | Cell Number | Medium (mL) | DNA (µg) | Transfection Reagent (µL) | DMEM (µL) |
|--------------|-------------------------|-----------------------|-------------|----------|---------------------------|-----------|
| 6 - Well | 10 | 9.5 x 10 ⁵ | 1.0 | 4 | 7.55 | 200 |
| 60 mm | 25 | 2.5 x 10 ⁶ | 3.0 | 7-12 | 15-24 | 300 |
| 100 mm | 55 | 5.5 x 10 ⁶ | 8.0 | 15-25 | 35-45 | 500 |
| 150 mm | 150 | 15 x 10 ⁶ | 20 | 45-75 | 60-80 | 1000 |

Protocol (example of producing virus by transfecting 293T cells in **one 100mm petri dish**.)

Refer to Table 1 for other plates or dishes)

Notice: Transfect cells at **70-90% confluency** for high transfection efficiency and low toxicity. Before plating, pipet the 293T cell suspension vigorously to achieve **single cells**. Plate the cells in **evenly** distribution manor.

Day 1: Plating cells

1. **Plating: 18 to 24 hours** prior to transfection, coat 100mm dishes with **6 mL 1x Gelatin** (101Bio, Cat. #: P910) for 60 min. Aspirate gelatin, and plate **5.5 x 10⁶** fast growing 293T cells.

Day 2: Transfection

2. **Change medium: 2 hours before** transfection, remove culture medium and add **8 mL fresh** complete culture medium (with 10% serum / **without antibiotics**)
3. In **tube 1** add:
7.5 µg Retrovirus vector DNA (Containing your gene of interest)
15 µg Retrovirus Packaging Mix (101Bio, Cat. # P905P)
500 µL DMEM (serum-free, High Glucose)
Pipet up and down to mix well
4. In **tube 2** Add:
45 µL 293T Transfection Reagent for Retrovirus Packaging
500 µL DMEM (serum-free, High Glucose)
Gently pipet up and down to mix well
5. Incubate at room temperature (**20–25°C**) for **3 min**.
6. Add **tube 2** into **tube 1**, pipet up and down several times. **Vortex** for **10 seconds**.
7. Incubate for **15 minutes** at **room temperature**.
8. Add the incubated mixture **drop-wise** to the cells, and **gently swirl** the plate to disperse evenly in the plate.
9. Return the cells to 37°C incubator with 5% CO₂.

Day 3: Add Titer-Up to increase virus titer

Add **20 µl** of Retrovirus **Titer-Up** (101Bio, Cat. #: 909) to the medium. Return the plates to the cell culture incubator.

Day 4: Collect virus

10. Collect virus supernatant twice at **48 and 72 hours** post transfection into a 50mL sterile FALCON tube. Centrifuge at **3,000rpm** for **15 minutes** at **4°C** to remove cell debris. Filter the clear supernatant through **0.45 µm** syringe tip filter.
11. The filtered clear supernatant is virus soup. Use it immediately or aliquot into sterile 1.5-mL tubes and store at **-80°C**, for **up to 3 months**.
12. (**Optional**) To concentrate virus, add **5X Retrovirus Concentration Solution** (101Bio, Cat. #: P905C) to the viral supernatant (volume of Retrovirus Concentration Solution vs. volume of viral supernatant = 1:4) and mix thoroughly. Put the mixture to 4°C refrigerator overnight and spin the virus pellet down next day. Please refer our P905C user manual for details.

-- The end --