


Detection of Alzheimer's Novel Biomarkers

pT-217, pT-181, pT-231 and Total Tau FAST ELISA Test

Biomarkers

Novel p-tau biomarkers: total tau, phosphorylated at threonine-181 (pT-181), threonine-217 (pT-217) and at threonine-231 (pT-231), increase more prominently in preclinical Alzheimer, and can accurately diagnose Alzheimer's Disease.

Locations of these three p-tau biomarkers




FLEXIBLE 80 uL of sample is enough for detection.

ACCURATE detect as low as 0.1 ng/mL.

SIMPLE combine samples and detections in one step.

TIME-SAVING use 2 hours instead of 5-6 hours

Our ELISA Standard Curves



Related ELISA kits

Fast Human Phospho-Tau Thr217 ELISA (T3293)

Fast Human Phospho-Tau Thr181 ELISA (T3294)

Fast Human Phospho-Tau Thr231 ELISA (T3296)

MAPT/Tau (total) ELISA (Cat.# T3295)

<https://www.101bio.com/ELISA.php>

Please find our manuals at www.101Bio.com