


Universal DNA Marker covers 100 bp - 10,000 bp

1 marker good for all DNA gel running



Cat.#: W2583-100 (100 rxn), W2583-500 (500 rxn)

List of our DNA Markers

Cat.#	DNA Marker	Size
W2583-100	Universal DNA Marker	500 µL (100 rxn)
W2583-500	Universal DNA Marker	2,500 µL (500 rxn)
W108-100	1 kb DNA Ladder 2.0	600 µL (100 rxn)
W108-500	1 kb DNA Ladder 2.0	3,000 µL (500 rxn)
W0636-100	100 bp DNA Ladder	500 µL (100 rxn)
W109-100	50/100bp DNA Ladder	600 µL (100 rxn)

