

## Retrovirus Packaging Service Requisition Form

### Customer information

|        |  |           |  |
|--------|--|-----------|--|
| Name   |  | Institute |  |
| E-mail |  | Telephone |  |

### Packaging service request

| Cat. # | Titer                        | Volume | Quantity Needed | Name of the virus |
|--------|------------------------------|--------|-----------------|-------------------|
| S502H  | > 1 x 10 <sup>8</sup> IFU/mL | 100 µL |                 |                   |
| S502HL | > 1 x 10 <sup>8</sup> IFU/mL | 200 µL |                 |                   |
| S502U  | > 1 x 10 <sup>9</sup> IFU/mL | 100 µL |                 |                   |
| S502UL | > 1 x 10 <sup>9</sup> IFU/mL | 200 µL |                 |                   |

### Vector information

|   | Vector name | Gene size (kb) | Volume (µL) | Concentration (µg/µL) |
|---|-------------|----------------|-------------|-----------------------|
| 1 |             |                |             |                       |
| 2 |             |                |             |                       |
| 3 |             |                |             |                       |
| 4 |             |                |             |                       |
| 5 |             |                |             |                       |

We require 20 µg endotoxin-free DNA for S302, and 40 µg endotoxin-free DNA for S302L and S302H. All custom virus production services should come with information on the plasmids which have to be packaged. Vector must be able to produce virus (no internal poly A signal, no toxic genes, no unusual 2nd structure, or 5' to 3' LTR size over 8kb.)

**Packaging system:** \_\_\_\_\_ **Packaging vectors provided?** Yes No

### Other requirement

Please send filled form by email ([support@101Bio.com](mailto:support@101Bio.com)) or fax (1-888-871-2390). Contact us to arrange plasmid shipment. We will response within 1 business day.